

CO ANALYSER
OF EXPIRED AIR



CATCH
YOUR BREATH

TABATABA^{v2}

**MEASUREMENT OF CO
IN EXPIRED AIR :**

- ▶ Help to stop smoking
- ▶ Intoxication level indicator light
- ▶ Instantaneous result in PPM,
% HbCo and % foetal HbCo

CE
0459



LCD screen
Instant result in PPM*,
% HbCo* and % foetal HbCo*

Disposable trapper
with no-return valve and
single-use mouthpiece



**SMOKING KILLS
...SO DOES
DOING NOTHING!**

Get into the habit of having your «**CO calculator**» to hand. The **TABATABA v2®** quickly and reliably tests levels of expired CO in active and passive smoking.

Indicator light
Green to red CO intoxication

*PPM parts per million
% HbCO: percentage of CO trapped in the haemoglobin
%HbCOF: percentage of CO trapped in the foetal haemoglobin

TABATABAv2®

THE +

- ▶ Easy to use
- ▶ Awareness for active and passive smoking cessation
- ▶ Accompaniment to stop smoking
- ▶ Helps motivation
- ▶ Low hydrogen, H₂S, SO₂, NO₂ interference
- ▶ Small-sized pipettes, trappers
- ▶ No technical knowledge required

YOUR NEW ALLY AGAINST ACTIVE AND PASSIVE SMOKING

Easy to use, accurate and reproducible, the **TABATABA v2®** is a measurement device intended for health professionals (tobaccologist, occupational health worker, associations, maternity departments, lung specialist, chemist, etc.) manufactured by FIM MEDICAL, leader in the screening field.

It measures the level of CO in expired air and provides an accompaniment for stop smoking programmes. Widely used by doctors, it is also a wonderful teaching instrument, its power of awareness to stop smoking smoking indisputable.

The **TABATABA v2®** is also useful during such public events as World No Tobacco day.

STEP 1



Breathe in

STEP 2



Hold your breath

STEP 3



Place the pipette
in your mouth and exhale slowly
and completely

STEP 4




Wait 15 seconds
and the result is displayed

YOUR SUPPLIER

37%
of smoking
adolescents aged
11-15 claim they are
addicted to nicotine.
2008 indicator
«les jeunes et le tabac» French
Cardiology Federation.

17,8%
of pregnant
women still
smoke during the
3rd trimester of
pregnancy.
(source INPES)

**1 cigarette
less**
equals 11 more
minutes of life.



ESSENTIAL ACCESSORIES

- ▶ Not sold with device:
 - Single-use pipettes (box of 250)
 - Single-use trappers (box of 250)
- ▶ 2 batteries included

**1
YEAR
WARRANTY**
TABATABA®



51 rue Antoine Primat - CS 60194
69625 Villeurbanne Cedex France

Tél +33 (0)4 72 34 89 89
Fax +33 (0)4 72 33 43 51
contact@fim-medical.com

CE
0459